



Recombinant *Drosophila virilis* Frizzled (fz)

Product Code	CSB-CF630520DMP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q24760
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila virilis</i> (Fruit fly)
Sequence	AQRYDQTPLDASSYYRSDLTGSSASSLDGFPHHNRCEPITISICKNIPYNMTIM PNLIGH TKQEEAGLEVHQFAPLVKIGCSADLQLFLCSLYVPVCTILERPIPPCRSLCESAR VCETL MKTYNFNWPENLECSKFPVHGGEDLCVAENTTASSSTPAPTRSAPKVTTTRKH QISVDSPH RNIGFVCPVQLKTPLGMGYELKVGGKDLHDCGAPCHAMFFPERERTVRLRYWV GSWAAICV ASCLFTVLTFLIDSSRFYPERAIVFLAVCYLVVGCAYVAGLGAGDSVSCREPF PPPVKL GRLQMMSTITQGHRQTTACTVLFMALYFCCMAAFWWSCLAFAWFLAAGLK WGHEAIENK SHLFHLVAWAVPALQTISVLALAKVEGDILSGVCFVGGQLDTHSLGGFLILPLCIY LSIGA LFLAGFISLFRIRTVMKTDGKRTDKLERLMLRIGFFSGLFILPALGLLGCLFYEY YNFD EWMIQWHRDICKPFSIPCPAARPPGTPEARPIFQIYMVKYLCMMLVGVTSVW LYSSKTM VSWRNFVERLQGKEPRTRAQAYV
Source	in vitro E.coli expression system
Target Names	fz
Protein Names	Recommended name: Frizzled Alternative name(s): Frizzled-1 Short name= dFz1
Expression Region	18-580
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.