



Recombinant *Drosophila melanogaster* Frizzled (fz)

Product Code	CSB-CF325586DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P18537
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	<p>YDQSPLDASPYRSGGGLMASSGTELDGLPHHNRCEPITISICKNIPYNMTIMP NLIGHT KQEEAGLEVHQFAPLVKIGCSDDLQLFLCSLYVPVCTILERPIPPICRSLCESAR VCEKLM KTYNFWPENLECSKFPVHGGEDLCAENTTSSASTAATPTRSVAKVTTRKH QTGVESPH RNIGFVCPVQLKTPLGMGYELKVGGKDLHDCGAPCHAMFFPERERTVLRVWV GSWAAVCV ASCLFTVLTFLIDSSRFYPERAIVFLAVCYLVVGCAYVAGLGAGDSVSCREPF PPPVKL GRLQMMSTITQGHRQTTCTVLFMALYFCCMAAFWWSCLAFAWFLAAGLK WGHEAIENK SHLFHLVAWAVPALQTISVLALAKVEGDILSGVCFVGQLDTHSLGAFILPLCIYL SIGA LFLAGFISLFRIRTVMKTDGKRTDKLERLMLRIGFFSGLFILPAVGLLGCLFYEY YNFD EWMIQWHRDICKPFSIPCPAARAPGSPEARPIFQIFMVKYLCSMLVGVTSVW LYSSKTM VSWRNFVERLQGKEPRTRAQAYV</p>
Source	in vitro <i>E.coli</i> expression system
Target Names	fz
Protein Names	Recommended name: Frizzled Alternative name(s): Frizzled-1 Short name=dFz1
Expression Region	19-581
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.