



Recombinant Arabidopsis thaliana MLO-like protein 10 (MLO10)

Product Code	CSB-CF863793DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9FKY5
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	<p>MATRCFWCWTTLLFCSQLLTGFARASSAGGAKEKGLSQTPTWAVALVCTFFIL VSVLLEK ALHRVATWLWEKHKNSLLEALEKIKAEMLILGFISLLLTFGGEQYILKICIPEKAAAS MLP CPAPSTHDQDKTHRRRLAAATTSSRCDEGHEPLIPATGLHQLHILLFFMAAFHI LYSFIT MMLGRLKIRGWKKWEQETCSHDYEFSDIPSRFRLTHETSFVRQHSSFWTKIPF FFYAGCF LQQFFRSVGRDYLTLRHGFIAAHLAPGRKFDQKYIKRSLEDDFKVVVVISPL LWASFV IFLLNVNGWEALFWASILPVLILAVSTKLQAILTRMALGITERHAVVQGIPLVHG SDK YFWFNRPQLLLHLLHFALFQNAFQLTYFFWVWYSFGLKSCFHTDFKLIVIKLSL GVGALI LCSYITLPLYALVTQMGSNMKKAVFDEQMAKALKKWHMTVKKKKGKARKPPT ETLGVSDT VSTSTSSFHASGATLLRSKTTGHSTASYMSNFEDQSMSDLEAEPLSPEPIEGH TLVRVGD QNTEIETGDISPGNQFSFVKNVPANDID</p>
Source	in vitro E.coli expression system
Target Names	MLO10
Protein Names	Recommended name: MLO-like protein 10 Short name= AtMlo10
Expression Region	1-569
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.