



Recombinant Mouse Frizzled-9 (Fzd9)

Product Code	CSB-CF870827MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9R216
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p>LEIGRFDPERGRGPAPCQAMEIPMCRGIGYNLTRMPNLLGHTSQGEAAAQLA EFSPLVQY GCHSHLRFFLCSLYAPMCTDQVSTPIACRPMCEQARLRCAPIMEQFNFGWP DSLDCARL PTRNDPHALCMEAPENATAGPTEPHKGLGMLPVAPRPARPPGDSAPGPGSG GTCDNPEKF QYVEKSRSCAPRCGPGVEVFWSSRRDKDFALVWMAVWSALCFFSTAFTVFTF LLEPHRFQY PERPIIFLSMCYNVYSLAFLIRAVAGAQSVACDQEAGALYVIQEGLENTGCTLVF LLLYY FGMASLWVVLTLTWFLAAGKKWGHEAIEAHGSYFHMAAWGLPALKTIVVL TLRKVAGD ELTGLCYVASMDPAALTGFVLVPLSCYLVLGTSFLLTGFVALFHIRKIMKTGGTN TEKLE KLMVKIGVFSILYTPATCVIVCYVYERLNMDFWRLRATEQPCTAATVPGGRR DCSLPGG SVPTVAVFMLKIFMSLVGITSGVWVWSSKTFQWQSLCYRKMAAGRARA CRTPGGYG RGTCHYKAPTIVLHMTKTDPSLENPTHL</p>
Source	in vitro E.coli expression system
Target Names	Fzd9
Protein Names	Recommended name: Frizzled-9 Short name= Fz-9 Short name= mFz3 Short name= mFz9 Alternative name(s): CD_antigen= CD349
Expression Region	24-592
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.