



Recombinant Human Frizzled-9 (FZD9)

Product Code	CSB-CF009112HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O00144
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	LEIGRFDPERGRGAAPCQAVEIPMCRGIGYNLTRMPNLLGHTSQGEAAAELAE FAPLVQY GCHSHLRFFLCSLYAPMCTDQVSTPIACRPMCEQARLRCAPIMEQFNFGWP DSLDCARL PTRNDPHALCMEAPENATAGPAEPHKGLGMLPVAPRPARPPGDLGPGAGGS GTCENPEKF QYVEKSRSCAPRCGPGVEVFWSRRDKDFALVWMAVWSALCFFSTAFTVLTFL LEPHRFQY PERPIIFLSMCYNVYSLAFLIRAVAGAQSVACDQEAGALYVIQEGLENTGCTLVF LLLYY FGMASSLWWWVLTLTWFLAAGKKWGHEAIEAHGSYFHMAAWGLPALKTIVILT LRKVAGD ELTGLCYVASTDAAALTGFVLVPLSGYLVLGSSFLLTGFVALFHIRKIMKTGGTN TEKLE KLMVKIGVFSILYTVPATCVIVCYVYERLNMDFWRLRATEQPAAAAGPGGRR DCSLPGG SVPTVAVFMLKIFMSLVVGITSGVWVWSSKTFQTWQSLCYRKIAAGRARAKAC RAPGSYG RGTHCHYKAPT VVLHMTKTDPSLENPTHL
Source	in vitro E.coli expression system
Target Names	FZD9
Protein Names	Recommended name: Frizzled-9 Short name= Fz-9 Short name= hFz9 Alternative name(s): FzE6 CD_antigen= CD349
Expression Region	23-591
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.