



Recombinant Arabidopsis thaliana MLO-like protein 4 (MLO4)

Product Code	CSB-CF522818DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O23693
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	<p>MEHMMKEGRSLAETPTYSVASVVTVLVFCFLVERAIYRFGKWLKKTRRKALF TSLEKMK EELMLLGLISLLLSQSARWISEICVNSSLFNSKFYICSEEDYGIHKKVLEHTSST NQSS LPHHGIHEASHQCQGHGREPFVSYEGLEQLLRFLFVLGITHVLYSGIAIGLAMSKI YSWRK WEAQAIIMAESDIHAKKTKVMKRQSTFVFHHASHPWSNNRFLIWMLCFLRQFR GSIRKSD YFALRLGFLTKHNLPTYNFHMVMVRTMEDEFHGIVGISWPLWVYAIVCICINV HGLNMY FWISFVPAILVMLVGTKLEHVVSKLALEVKEQQTGTSTNGAQVKPRDGLFWFGK PEILLRL IQFIIFQNAFEMATFIWFLWGIKERSCFMKNHVMISSRLISGVLVQFWCSYGTVP LNVIV TQMGSRHKKAVIAESVRDSLHSWCKRVKERSKHTRSVCSLDTATIDERDEMT VGTLSRSS SMTSLNQITINSIDQAESIFGAAASSSSPQDGYTSRVEEYLSETYNNIGSIPPLN DEIEI EIEGEEDNGGRGSGSDENNGDAGETLLELFRRT</p>
Source	in vitro E.coli expression system
Target Names	MLO4
Protein Names	Recommended name: MLO-like protein 4 Short name= AtMlo4
Expression Region	1-573
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.