



Recombinant Rabbit Alpha-1D adrenergic receptor (ADRA1D)

Product Code	CSB-CF001387RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O02666
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	MTFRDLLSVTFEGPRPDISAGGSGAGGGAGAGAGAGDTASSESPAVGGVPG AAGGGGGGS VVGAGSGEDNRSSAGEPGGAGGGGGEVNGTAAVGGGLVVSQAQSVGVGVFLAA FILTAVAGNL LVILSVACNRHLQTVTNYFIVNLAVADLLLSATVLPFSATMEVLGFWAFGRAFC DVWAAV DVLCCASILSLCTISVDRYVGVVRHSLKYPAIMTERKAAAILALLWAVALVSMG PLLGW KEPVPPDERFCGITEEVGYAVFSSLCSFYLPMAVIVVMYCRVYVVARSTTRSLE AGVKRE RGKASEVVLRIHCRGAASGADGAPGTRGAKGHTFRSSLSVRLKFSREKAA KTLAIVVG VFVLCWFPPFFVLPLGSLFPQLKPSEG VFKVIFWLGYFNVCNPLIYPCSSREF KRAFLR LLRCQCRRRRRRPLWRVYGHWRASAGGGPHPCALSAGAALPGAALALT AAPASSAA APEGQAAGAGRRKPPCAFREWRLGPLRRPTTQLRAKVSSLSHKIRAGGAQR AEAACALR SEVEAVALSVARDVAEDNTCQAYELADYRNLRETDI
Source	in vitro E.coli expression system
Target Names	ADRA1D
Protein Names	Recommended name: Alpha-1D adrenergic receptor Alternative name(s): Alpha-1D adrenoreceptor Short name= Alpha-1D adrenoceptor
Expression Region	1-576
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.