



Recombinant Human XK-related protein 7 (XKR7)

Product Code	CSB-CF685936HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5GH72
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAAKSDGAAASASPDPEGAAGGARGSAGGRGEAAAAAGPPGVVAGGGPGP RYELRDCCWV LCALLVFFSDGATDLWLAASYLQHQHTYFSLTLLFVLLPSLVVQLLSFRWFVY DYSEPA GSPGPAVSTKDSVAGGAAISTKDSAGAFRTKEGSPEPGPQPAPSSASAYRRR CCRLCIWL LQTLVHLLQLGQVWRYLRALYLGLQSRWRGERLRRHFYWQMLFESADVSM LLETFLLRS APQLVLQLSLLVHRGGAPDLLPALSTSASLVSLAWTLASYQKVLDRSRDDKRP LSYKGAV AQVLWHLFSIAARGLAFALFASVYKLYFGIFIVAHWCVMTFWVIQGETDFC MSK WEEIY NMVVGIIYFCWFNVKEGRSRRRMTLYHCIVLLENAALTGFWYSSRN FSTDFYS LIMVCV VASSFALGIFFMCVYYCLLHPNGPMLGPQAPGCIFRKASEPCGPPADAI TSP RSLPRTT GAERDGASAGERAGTPTPPVFQVRPGLPPTPVARTLRTEGPVIRIDL PRKKYP AWDAHFI DRRLRKTILALEYSSPATPRLQYRSVGTSGELLE YETTV
Source	in vitro E.coli expression system
Target Names	XKR7
Protein Names	Recommended name: XK-related protein 7
Expression Region	1-579
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.