



# Recombinant *Xenopus laevis* Hyaluronan synthase-related protein (has-rs)

<b>Product Code</b>	CSB-CF524498XBE
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O57428
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Xenopus laevis</i> (African clawed frog)
<b>Sequence</b>	<p>MENTTDPENIPVSKPKYPTIRRILSQTFRILLIFSITTAYVLGYQALCHQGGLLITFG  LYG  AAMLLHLLMQGIFANLEIRRIEKRDGVCSFKKTVALTITGYQENPDYLRQCLES  CKGMKY  PKDKLKIILVIDGNNEEDVYMMEIFKEVFHGEDVGTYYVWQENYHTWNIPSEESE  DSSSEI  SSFPWKNEGIQMVEELVRTKRCVCIMQQWGGKREVMYTAFRALGTSVDFILV  CNSDIKLD  KMATVELVKVLEDDDKNGAVGGDVRVWNRHDSFISFMSSLRYWMVFNMEIA  CQSYFDSVT  YIRGSLGMYRNDILQAFLEFWYNKTFLGTRCPIGDDRFLTNRVLSMGYRTKYS  HKSCAYA  PCQYLRWLNQQTPWARSYFRMWFCNAQWWHQHHIWMTYESATGIFFPFV  TAVLIRLMYS  SSLCNIVWLFLCIQIMSLLSLYASWQSKKLSMVLMSLYSTLYIIWLLPCQLVALL  TIAK  SDWGTSGRKKVVNNYVPLFSLSIWAAVLLGGLCYSMYIGCRKDWSKPQANRE  LYHLLYGC AGYMAYWVLMTVIYCVSGSCCKMRSQAVPQTHDITSLSVSLLV</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	has-rs
<b>Protein Names</b>	Recommended name: Hyaluronan synthase-related protein Alternative name(s): Hyaluronan synthase-related sequence Short name= xHAS-rs
<b>Expression Region</b>	1-583
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.