



Recombinant *Arabidopsis thaliana* Uncharacterized membrane protein At3g27390 (At3g27390)

Product Code	CSB-CF812740DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8GUM4
Product Type	Transmembrane Protein
Immunogen Species	<i>Arabidopsis thaliana</i> (Mouse-ear cress)
Sequence	MEPPIGFRASLFQFLLFLPYFIGLLFLGFIKIVLCPLVCLVVTIGNSAVILSLLPVH IV WTFYSIVSAKQVGPILKIFLCLCLPAAIILWPIVGILGSVLGGALYGFFSPIFATFD AVG EGKPYQFFHCFYDGTWSTMQRSFTVVRDFKDVCFHSYFSLMDELKQSCPDR KYYEIRLLQ LPGALVVSVLGILVDPPVISLVAICKSPYMLFKGWHRLFHDLIGREGPFLETMCV PIAGL AILLWPLAVTGAVIGSVISSIFLGAYAGVVSQESSFYGLCYIVASVSIYDEYST DILD LPEGSCFPRPKYRRKDEEPTPFSGVPVRLGSVKNASSMRGGSVRVPMIDIKPL DLLNELF VECRRYGEVLATKGLINSKDIEEARSSKGSQVISVGLPAYGLLYEILRSVKANSS GLLLS DGVTEITTMNRPKDVFFDWFLNPFLILKEQMKATNLSEEEEEYLGRVLVLLFGDP ERLKSS NAISAPPLTERKRAELDAFARMQGLTKTVSRYPTFRRHFVALVKKLSEDLDL KDNNSA KDESITEPPAPVKIISRIFSQRSFRRKGSVNGSDQESQKGVSRNVDIV
Source	in vitro E.coli expression system
Target Names	At3g27390
Protein Names	Recommended name: Uncharacterized membrane protein At3g27390
Expression Region	1-588
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.