



Recombinant Gorilla gorilla gorilla Muscarinic acetylcholine receptor M3 (CHRM3)

Product Code	CSB-CF865050GGZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9N2A3
Product Type	Transmembrane Protein
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Sequence	<p>MTLHNNSTTSPLFPNIISSSWIHSPSDAGLPPGTDTHFGSYNVSRAAGNFSSPD GTTDDPL GGHTVWQVVFI AFLTGILALVTIIGNILVIVSFKVNKQLKTVNNYFLLSLACADLIIG VI SMNLFTTYIIMNRWALGNLACDLWLAIIDYVASNASVMNLLVISFDRIYFSITRPLT YRAKR TTKRAGVMIGLAWVISFVLWAPAILFWQYFVGKRTVPPGECFIQFLSEPTITFGT AIAAF YMPVTIMTILYWRIYKETEKRTKELAGLQASGTEAETENFVHPTGSSRSCSSYE LQQQSM KRSNRRKYGRCHF WFTTKSWKPSSEQMDQDHSSSDSWNNNDAAASLENSA SSDEEDIGSE TRAIYSIVLKLPGHSTILNSTKLPSSDNLQVPEEELGMVDLERKADKLQAQKSV DDGGSF PKSFSKLPIQLES AVDTAKTSDVNSSVGKSTATLPLSFKEATLAKRFALKTRSQI TKRKR MSLVKEKKAQTL SAILLAFIITWTPYNIMVLVNTFCDSCIPKTFWNLGYWLCYI NSTVN PVCYALCNKTFRTTFKMLLLCQCGKKRRKQQYQQRQSVIFHKRAPEQAL</p>
Source	in vitro E.coli expression system
Target Names	CHRM3
Protein Names	Recommended name: Muscarinic acetylcholine receptor M3
Expression Region	1-590
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.