



Recombinant Dictyostelium discoideum Frizzled and smoothened-like protein H (fslH)

Product Code	CSB-CF690806DKK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q556C6
Product Type	Transmembrane Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Sequence	<p>QRYLPVEGGKCEKYIGDGDGSGKRICNGYLANPDSVYVHNKTQQETLKDLRSLI NLELNN PSKECINPSNYKIMCAMMFPECIEINGTNIVKPITLPIYTCNSFCKEALVTCVNP TIAS CDGGTNLPIQLPYTPIEWVKYNLTIYGGVDDYRVNCTDPTLISDSGSSSEIEVG CVEPLI KRPTNDTKGDLEKGYFYVNSQCVINCPVTGMHPKSVWNQIFKINDVLSSISLAC TLILLF TFGILNPKLNRFDKKNLFFIAGVFGMSVSVGLIAANGSEKTVCPTEPERYAVNTD RVCVAS GFLVHFSALFAILWWTIGLADVYYGIKIVGKKIKIKVRYLLATLTISLAFTLVPLG TGQ YQAGLSNVMCFLKDEIYQSMTFFVPLGICLTMGTILMILVMREIYVIVKSNSTSS SFSSS SSKSKSKSKSSDSISYLKLVKPKMLNIIIFYFTFLYFLFVRVINSRYQEYEDSAI PYML CLAKGGGDSCRLKGPSAGSLGYFAYCLRIYGIYLFIIISFLSSRTIKIWKESIILNNA FVT PIIKFIDSSFSNRFSSSKNTSTTQNSTLNNTESDTSKRGNSSAVSINLESRYNT DDDDL</p>
Source	in vitro E.coli expression system
Target Names	fslH
Protein Names	Recommended name: Frizzled and smoothened-like protein H
Expression Region	22-621
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.