



Recombinant Rat Choline transporter-like protein 3 (Slc44a3)

Product Code	CSB-CF731823RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6AY92
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>MGYSVVAGAAGRLLFGYDSFGNVCGKRNSPVEGAPLSGQDMTLKKHVFFMN ACNLEVKDR GRGPMALCVSSCPEKQLDTLEEVQLFANINGSFLCVYSLNSFNQTQSPSADRL CPRLPVP PSKPFPLFNRCIPQTPECYSLFASVLLNDADALHRILSGIMAGRDTILGLCAFVF ALSLA MLFTFRFISTLLVHIIISLVILGLLFVCGVFWWLYDYTNDLSTELDTEKENMKCM LAFA VITTVTVVLLALIFTLRKRVKLTVELLRVTNKAISRCPFLLLQPLWTFAILVFFWV LWV AVLLSLGTAGTAQVMEGGQVEYKPLSGIRYLWWYHLIGLIWTSEFILTQQRMT VAGAMVT CYFNRNQNDPPARPILSSLSTLFCYHQGTAVKGSLLLVTRIPRVIFMYIYSTVK ERHSA WPRGEFRCSYCGLWCLTNYPYHLNQDAYAAAINGTDFCTSAKDAHTIISKNS SHLTSVN CFGNFIIFLGKVLVVCFSVFGGLMAFNYNRALQVWAIPLLLVAFFAYLAHSFLS VFETV LDTLFLCFAVDLETNDGSSEKPYFMDPGFLSFVKRTDHFNNARSQGHKNSLP NEEGTELR PIVR</p>
Source	in vitro E.coli expression system
Target Names	Slc44a3
Protein Names	Recommended name: Choline transporter-like protein 3 Alternative name(s): Solute carrier family 44 member 3
Expression Region	1-604
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.