



Recombinant Dictyostelium discoideum Frizzled and smoothened-like protein M (fsIM-1)

Product Code	CSB-CF684822DKK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q556N7
Product Type	Transmembrane Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Sequence	<p>QTIYPIDPSGKCEQYIGDTPSSPCSKFLNNLDSIYVSSNATQQNAMKKLDEYFG ILGAIG TSGCKSDSLTYKTLCSIYLPGESFTNNQTNITIAIPKRICYNTCNNVTTKCGVP KFYFS CDQIEPVSGLPMPFLNYSEFNLTNYDVGPNPYTVQCYGPLSDNTMVNLDSSY CPPPLFYH NSTDHDADYDRGYLFVSETSNCVVPNPVPLYTNEQWDQLYNLSNTLAVLSTF GSLYLLVT FIILNPKVTSFDRMYGFFNGSVFMMSLSGVILFIAGGPRALIKDGGARISVFEDP LCSST GFIFQLFAINAILFWAYMGFDLWWRVKYITKPLNIQKYYVPIAFTISFIFSVIPLAT KNY RMVRGNIHCWVHKAVLQNTLFFGPLGLTLTISTGFIGLVIYEIYKIVKATGRGGIM KLEI KPILNIVLIYFSFVYIFAFNFHNDNNSKNTYGSIDEFFQCTLESDDPSKCTVGGP SIGSL GYFIYCIRIYGIYCFFLQGLNERAFKIWKRSIVFNNRFILFIKVKLFSMDNNSPSES GNS STTAGTSTTINNSNINKKNNNSKPTLSTMDSNFAFSKNNDSDFSDFDDYDPYHKK QNDIEIG SVNIK</p>
Source	in vitro E.coli expression system
Target Names	fsIM-1
Protein Names	Recommended name: Frizzled and smoothened-like protein M
Expression Region	19-623
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.