



# Recombinant Bacteroides thetaiotaomicron Membrane protein insertase YidC (yidC)

<b>Product Code</b>	CSB-CF805825BDU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q8AA76
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Bacteroides thetaiotaomicron (strain ATCC 29148 / DSM 2079 / NCTC 10582 / E50 / VPI-5482)
<b>Sequence</b>	<p>MDKNTITGLVLIGILLVGFSYLSRPSEEQIAAQKKYYSIAVVQQQQEALKAKTE  AALAN  ENKGAAAAADSSALFFNAMHGTDISKVSIQNSVAEITFTTKGGRVYSAMLKEYK  GQDKTNP  VVLFDGDDATMSFNFYNKQGAIQTKDYIFEAVNKTSSVTMRLAADNASYIDFI  YTLKPN  SYLMNFEIKATGMEGKLASTEYVDIDWTQRARQLEKGFTYENRSELTYKVKG  DNVDNLS  AAKDDEKDLGNTAIDWVAFKNQFFSSVFIADQDFNKVSVKSRMEQQGSGYIKD  YSAEMST  FFDPSGKQPTEMYFYFGPNHFKTLKALDKGRTEKWELNRLVYLGWPLIRWIN  QFITINVF  DWLSGWGLSMGIVLLILTIMVKVVVYPATWKTYMSSAKMRVLKPKIDEINKKYP  KQEDAM  KKQQEVMSLYSQYGVSPMGCLPMLLQFPILMALFMFVPSAIELRQQSFLWA  DDLSTYDA  FITFPFHIPFLGNHLSLFCLLMTVTNINLNTKYTMTMQDTGAQPQMAAMKWMY  LMPIMFL  FVLNDYPSGLNYYYFVSTLISVGTMILLRKTDETKLLAILEAKKKDPKQMKKTG  FAARL EAMQKQQEQLQQQKQNKR</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	yidC
<b>Protein Names</b>	Recommended name: Membrane protein insertase YidC Alternative name(s): Foldase YidC Membrane integrase YidC Membrane protein YidC
<b>Expression Region</b>	1-618
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.