



Recombinant *Saccharomyces cerevisiae* Protein DFG16 (DFG16)

Product Code	CSB-CF859499SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q99234
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MIIRLHFYLLTLVYHLGLVGAYEKAARKRIQPPDLIPGPPGHKLGDERPPHYD HRPPYK KHIDNIPAYNLTDLIDDKLLNKYENSCTVNVLTGGFISLASNSWHLRAYNYTLNY PSFLI RCDNGSANPNFQSHVLQDFVYDINNKFNVQDDSSKYIGKDPFPLGMIMITFASG CICVATW MLFLVLLLLPSDNHNRRNKVVHVYVLFSAIIRTVFLNETIAVIFDSQYHDDYQDA SQFES FIVETAPYKICELVANILSDINWIYIVHYLQSNYGKPTWNWIPFKMKKGTHIIITVG CFL SLADNILFANLLWRKNLVVLKVIFYKLIELLIYTIFISIICYFTWHNFAYILLPKTAEIN T DGKCKTKLRILWENYHETIPLLAYNILIFILFYFTTIFFAAFTKHVRGWTFFNVHLL KVL ITVNVWGLIGVLEKRELHISKKTVLGRKINNRRDKFFANPTVNYYGEDLGKHL SAITLNRD LNTTKSNTTSHDSSSLVGSPSPTWKSPIERIRDRRRRHKIMKSENKFGQNPSF GSKSNGK PNTKTTLSKYRQLLRKPRRKTSYEPKNGIGQNKEGSTVRPGADKHIRDSNYL ATDISDN ESMETELRTNHIYNYNSD
Source	in vitro E.coli expression system
Target Names	DFG16
Protein Names	Recommended name: Protein DFG16 Alternative name(s): Extracellular mutant protein 41 Zinc-regulated gene 11 protein
Expression Region	1-619
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.