



Recombinant Mouse Probable G-protein coupled receptor 153 (Gpr153)

Product Code	CSB-CF809076MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8K0Z9
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MSDERRLPSSAVGWLACGGLSLLANAWGILSVGAKQKKWKPLEFLLCTLAAT HMLNVAVP IATYAVVQLRRQRDPDYEWNEGLCKVFVSTFYTLTLATCFSVTSISYHRMWMVR WPVNYRL SNAKKQAVHTVMGIWMVSFILSALPAVGWHDTSERFYTHGCRFIVAEIGLGFG VCFLLLV GGSVAMGMVCTAIALFQTLATQVGHRADRRFTTVPTIVVEDAQQKRRSSIDGS EPARTSL QITGLVATIVVIYDCLMGFPVLVVSFSSLRADASAPWMALCVLWCSVTQALLP LFLWTC DRYRADLKAVWEKCVALMANDESDNETSLEGSISPDMVLERSLDYSYGGDF VALDRMAK YELSALEGLPQLYPLRPLQEDRMQYLQGAGRHRCGFPGGQPCFSPAGCSQ VPPTRRFSH DDADVWAAVPLPTFLPRWSSGEDLAALAHMLPAGSDRRRGSLLAFAEDAPP FRPRRRA ESLLSLQPSSLDGGPRHAQDSPPGSPRRRPGPGARSASVSLLPDAFALTAFE REPQALRR VPAPAQPFFPAARDSAEPAEVPTPPGGRTQRSQGRRAARTHVGPLQSSLSAS WGEPGGLHA AGCGSISSFLSSPSESSGYVTLHSDSLGSAS
Source	in vitro E.coli expression system
Target Names	Gpr153
Protein Names	Recommended name: Probable G-protein coupled receptor 153 Alternative name(s): G-protein coupled receptor PGR1
Expression Region	1-631
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.