



# Recombinant *Yarrowia lipolytica* pH-response regulator protein palH/RIM21 (RIM21)

<b>Product Code</b>	CSB-CF892268YAP
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q9UVF6
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) ( <i>Candida lipolytica</i> )
<b>Sequence</b>	MHSDAPTPVANELSEVSHLAEGDNASFSAGNFTVLQLPSDPSHCIDYAIPAGT LVIVDPN NPDNNRTVKLQAPAVFRPQCALGGTRAPAPRPEESFPDWSEYQKYHHNDRR DPFYGSVTP IAYTIAASTVTAWMLLILFLSRKPSPLFQKIAVLITAVSLTVFLAQATDTLESQYN EGY QNAYELRHKIMGGWAFRILQVITCVITWLARLQVVIRLFDVVPKINTRLAVVGSTLI FTNA TIWACLNLIPPWSQYVRNAKSVLPVFGALCSLLLEVFYLVVVVIYSISKRKYAYS RTSIV MAAISWLAMILPMVFIVFDIAHYWIAGWSDFIRWTADAAASVVVWEWTNVIVYQ ERREQR QSVLGRQVYRDEILDFKGDNGGGTVGGGRTKYPSRMEEDDVPRSSPNDHH FTSNIPTSA GEGQSFQFFKRARLPMYSRKIWKIARGESAASNTHYEHAIIEEEEEESIERNR RTPTVQ ENGEEDDEETYDEENDQYSQDNHSSVHSFESSRPSQHPVPSSGGTRAVGHT HFPLPGQSE GAHTTSPAAAAAAAEPEPEPVAGPSGGAAAHGDSDDDSDNSDDSSLASFTVI QQTGFSVD NQGVPEYDADSAPPTFEPIPGFHRQDYSDAKG
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	RIM21
<b>Protein Names</b>	Recommended name: pH-response regulator protein palH/RIM21
<b>Expression Region</b>	1-632
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.