



Recombinant Chicken Muscarinic acetylcholine receptor M3 (CHRM3)

Product Code	CSB-CF005383CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49578
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	<p>MLTHYQLCFQKRSSQNYTVPDPTSCFDVPPWTILCQRATMIMHNNSSALPLFS NVSSFWK RDSHGPGLLDEGASLIGSYDSPETTESFPFSTVETTNSLNLATIKDPLGGHAV WQVVLIA FLTGIIALVTIIGNILVIVSFKVKNQLKTVNNYFLLSLACADLIIGVISMNLFTTYIIMG HWALGNLACDLWLSIDYVASNASVMNLLVISFDYRFSITRPLTYRAKRTTKRAG VMIGLA WIISFVLWAPAILFWQYFVGKRTVPLDECIFIQLSEPIITFGTAIAAFYLPVTIMSIL YW RIYKETEKRTKELAGLQASGSEAETARFVHQTGSSRSLSSYELQRQSTKRSSR RKYRRCH FWLTMKSWEPNTDQGDQEHSSSDSWNNNDAAASLENSASSDEEDITAETRAI YSIVLKLP GHSAILNSTKLPSSDLNESADELQKSDTDSQEKKPKKLQPPKSIQDGGSFQK SFSKLPI QPGSAETATASDGISSVTKTSALPLSFKEATLAKKFALKTRSQITKRKRMSLIK EKCAA QTLSAILFAFIITWTPYNIMVLVNTFCDCVPKTVWNLGWLCYINSTVNPVCYAL CNKMF RNTFKMLLLCQCDKRKRKRRKQQYQQRQSVIFHKRIPREAS</p>
Source	in vitro E.coli expression system
Target Names	CHRM3
Protein Names	Recommended name: Muscarinic acetylcholine receptor M3
Expression Region	1-639
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.