



Recombinant *Xenopus laevis* Frizzled-3 (fzd3)

Product Code	CSB-CF009106XBE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O42579
Product Type	Transmembrane Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Sequence	<p>QKSMGHS�FACEPITLRMCQDLPNSTFMPNLLNHYDQQTAALAMEPFHPMV NLECSRDL RPFLCALYTPVCM EYGRMTLPCRKLCQRAYNECFKLMEMFGVPWPEEM ECS RFPDCDEPY PRIVDISLSGEPSEETPLAVQRDYGFWCPRELKIDPDLRSSF LGVRDCSPPCP HMYFRRE ELSFARYFIGVISIVCLSATLFTFLTLIDVTRFRYPERPIIFYAVCYMMVSLIFFIG FL LEDKVACNGANPSQYKASTVTQGSHNKACTMLFMVLYFFT MAGSVWWWILT I TWFLAAVP KWGSEAIEKKALLFHASAWGIPGTLTIILLAMNKIEGDNISGVCFVGLYDVHALR YFVLA PLCLDVVVGVSLLLAGIISLNRVRIEIPLEKENQDKLVKFMIRIGVFSILYLVPLL VV IG CYFYEQAYRGVWETT WVQERCREYHIPCPYKVTQT SRPDLILFLMKYLMLLVV GIPSVFW VGSKKTCFEWASFFHGRKKKAGVNESRQVLQEPDFAQSLLRDPNTPIVRKSR GTSTQGTS THASSTQLAMLDDQRSKAGSVQSKVSSYHGSLHRSRDGRYTPCSYRGIEERL PHGSMSHL TDHSRHSSTHRLNEQSHQGSIRDLSNPLAHISHGTS MNRVIEADATSA</p>
Source	in vitro E.coli expression system
Target Names	fzd3
Protein Names	Recommended name: Frizzled-3 Short name= Fz-3 Short name= Xfz3
Expression Region	17-664
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.