



Recombinant Rat Frizzled-8 (Fzd8)

Product Code	CSB-CF677711RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q498S8
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>ASAKELACQEITVPLCKGIGYNYTYMPNQFNHDTQDEAGLEVHQFWPLVEIQCS SPDLKFF LCSMYTPICLEDYKKPLPPCRSVCERAKAGCAPLMRQYGFAPDRMRCDRLP EQGNPDTL CMDYNRTDLTTAAPSPRRLLPPPQPGEQPPSGSGHSRPPGARPPHRGGSSR GSGDTAAAP PSRGGKARPPGGGAAPCEPGCQCRAPMVSSEHPLYNRVKTGQIANCAL PCHNPFFSQ DERAFTVFWIGLWSVLCFVSTFATVSTFLIDMERFKYPERPIIFLSACYLFSVSVG YLVRL VAGHEKVACSGGAPGAGGAGGAGGAATAGAGAAGAGASSPGARGEYEELG AVEQHVRVYET TGPALCTVVFLVYFFGMASSIWWVILSLTWFLAAGMKWGNEAIAGYSQYFHL AAWLVPS VKSIAVLALSSVDGDPVAGICYVGNQSLDNLRGFVLAPLVIYLFITMFLLAGFV SLFRI RSVIKQQGGPTKTHKLEKLMIRLGLFTVLYTPAAVVVACLFIYEQHNRPRWEA THNCPCL RDLQPDQARRPDYAVFMLKYFMCLVVGITSGVWVWSGKTLESWRALCTRCC WASKGAAVG AGAGGSGPGGGGPGGGGGHGGGGGSLYSDVSTGLTWRSGTASSVSYP KQMPLSQV</p>
Source	in vitro E.coli expression system
Target Names	Fzd8
Protein Names	Recommended name: Frizzled-8 Short name= Fz-8
Expression Region	28-684
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.