



Recombinant Mouse Frizzled-8 (Fzd8)

Product Code	CSB-CF726733MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q61091
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p> ASAKELACQEITVPLCKGIGYNYTYMPNQFNHDTQDEAGLEVHQFWPLVEIQCS SPDLKFF LCSMYTPICLEDYKKPLPPCRSVCERAKAGCAPLMRQYGFAPDRMRCDRLP EQGNPDTL CMDYNRTDLTTAAPSPRRRLPPPPPPGEGPPSGSGHSRPPGARPPHRGGSS RGSDDAAAA PPSRGGKARPPGGGAAPCEPGCQCRAPMVS SVSSERHPLYNRVKTGQIANCA LPCHNPFFS QDERAFTVFWIGLWSVLCFVSTFATVSTFLIDMERFKYPERPIIFLSACYLFVSV GYLVR LVAGHEKVACSGGAPGAGGAGGAGGAAAAGAGAAGAGASSPGARGEYEEL GAVEQHVRYE TTGPALCTVVFLLVYFFGMASSIWWVILSLTWFLAAGMKWGNEAIAGYSQYFH LAAWLVP SVKSI AVLALSSVDGDPVAGICYVGNQSLDNLRGFVLAPLVIYLFIGTMFLLAGF VSLFR IRSVIKQQGGPTKTHKLEKLMIRLGLFTVLYTPAAVVVACL FYE QHNRPRWEA THNCPC LRDLPDQARRPDYAVFMLKYFMCLVVGITSGVWVWSGKTLESWRALCTRC CWASKGAAV GAGAGGSGPGGSGPGGGGGHGGGGGSLYSDVSTGLTWRSGTASSVSY PKQMPLSQV </p>
Source	in vitro E.coli expression system
Target Names	Fzd8
Protein Names	Recommended name: Frizzled-8 Short name= Fz-8 Short name= mFz8
Expression Region	28-685
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Protein Length	Full Length of Mature Protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>