



Recombinant Mouse Frizzled-8 (Fzd8)

Product Code	CSB-CF726733MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q61091
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	ASAKELACQEITVPLCKGIGYNYTYMPNQFNHDTQDEAGLEVHQFWPLVEIQC SPDLKFF LCSMYTPICLEDYKKPLPPCRSVCERAKAGCAPLMRQYGFAPDRMRCRLP EQGNPDTL CMDYNRTDLTTAAPSPRRRLPPPPPPGEQPPSGSGHSRPPGARPPHRGGSS RGSMDAAAA PPSRGGKARPPGGGAAPCEPGCQCRAPMVSVSSERHPLYNRVKTGQIANCA LPCHNPFFS QDERAFTVFWIGLWSVLCFVSTFATVSTFLIDMERFKYPERPIIFLSACYLFVSV GYLVR LVAGHEKVACSGGAPGAGGAGGAGGAAAAGAGAAGAGASSPGARGEYEEL GAVEQHVRYE TTGPALCTVVFLLVYFFGMASSIWWVILSLTWFLAAGMKWGNEAIAGYSQYFH LAAWLVP SVKSI AVLALSSVDGDPVAGICYVGNQSLDNLRGFVLAPLVIYLFIGTMFLLAGF VSLFR IRSVIKQQGGPTKTHKLEKLMIRLGLFTVLYTVPAAVVACLIFYEQHNRPRWEA THNCPC LRDLPDQARRPDYAVFMLKYFMCLVVGITSGVWVWSGKTLESWRALCTRC CWASKGAAV GAGAGGSGPGGSGPGGGGGHGGGGGSLYSDVSTGLTWRSGTASSVSY PKQMPLSQV
Source	in vitro E.coli expression system
Target Names	Fzd8
Protein Names	Recommended name: Frizzled-8 Short name= Fz-8 Short name= mFz8
Expression Region	28-685
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.