



Recombinant Mouse Adhesion G-protein coupled receptor G1 (Adgrg1)

Product Code	CSB-CF812970MO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8K209
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Transmembrane Proteins
Immunogen Species	Mus musculus (Mouse)
Sequence	GSPREDFRFCGQRNQTQQSTLHYDQSSEPHIFVWNTTEETLTIRAPFLAAPDIP RFFPEPR GLYHFCLYWSRHTGRLHLRYGKHDYLLSSQASRLLCFQKQEQLKQGAPLIAT SVSSWQI PQNTSLPGAPSFIFSFHNAPHKVSHNASVDMCDLKKELQQLSRYLQHPQKAA KRPTAAFI SQQLSLESKLTSSVFLGDTLSFEEDRVNATVWKLPPPTAGLEDLHIHSQKEEE QSEVQAY SLLLPRAVFQQTRGRRRDDAKRLLVVDVDFSSQALFQDKNSSQVLGEKVLGIVVQ NTKVTNL SDPVVLTFFQHQPQPKNVTLQCWFVVEDPASSSTGSWSSAGCETVSRDTQTS CLCNHLTYF AVLMVSSSTEVEATHKHLYLTLLSYVGCVISALACVFTIAAYLCSRRKSRDYTIKVVH MNLIS AVFLLDVSFLLSEPVALTGSEAACRTSAMFLHFSLLAACLSWGMGLEGNLYRLV VEVFGTY VPGYLLKLSIVGWGFPVFLVTLVALVDVNNGPIILAVRRTPERVTYPSMCWIR DSLVSY VTNLGLFSLVFLFNLAAMLATMVVQILRLRPHSQNWPHVLTLLGLSLVLGLPWAL VFFSFA SGTFQLVILYLFSIITSFQGFLIFLWYWSMRFQAQGGPSPLKNNSDSAKLPISG STSSS RI
Source	in vitro E.coli expression system
Target Names	Adgrg1
Protein Names	Recommended name: G-protein coupled receptor 56 Alternative name(s): Serpentine receptor cyt28
Expression Region	26-687
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

**Tag Info**

N-terminal 10xHis-tagged

Protein Length

Full Length of Mature Protein

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.