



# Recombinant Donkey Follicle-stimulating hormone receptor (FSHR)

<b>Product Code</b>	CSB-CF842510DK
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q95179
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Equus asinus (Donkey) (Equus africanus asinus)
<b>Sequence</b>	<p>CHHQVCHYSNRVFLCQESKVTEIPSDLPRNALELRFVLTKLRVIPKGAFSGFGD  LKKIEI  SQNDVLEVIEANVFSNLPKLHEIRIEKANLLYIDHDAFQNLPNLQYLLISNTGIK  HLP  VHKIQSLQKVLLDIQDNINIHIVERNNSFMGLSFESMILRLSKNGIQEIHNCANFNGT  QLDE  LNLSDNNLEELPNDVVFQGASGPVILDISGTRIHSLPNYGLLENLKKLRARSTYN  LKKLPS  LEKFVALMEASLTYPHCCAFANWRQQTSELQTTCNKSILRQEVDMTQARGE  RVSLAEDD  ESMMYSEFDYDLCNEVVDVTCSPKPAFNPCEIMGYDILRVLWIFISILAITGN  IIVLV  ILITSQYKLTVPRFLMCNLAFAFADLCIGIYLLLIASVDIHTKSQYHNYAIDWQTGAG  CDAA  GFFTVEGSELSVYTLTITLERWHTITHAMQLECKVQLRHAASVMLVGVWIFGFG  VGLLPI  FGISTYMKVSIKLPMDIDSPSLQYVMSLLVLNVLAFVVICGCYTHIYLTVRNPNI  VSSS  SDTKIAKRMGILIFTDFLCMAPISFFGISASLKVALITVSKSKILLVLFYPINSCANP  FL  YAIFTKNFRRDFILLSKFGCYEMQAQTYRTETSSTGHISHPKNGPCPPTPRVT  NGANCT LVPLSHLAQN</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	FSHR
<b>Protein Names</b>	Recommended name: Follicle-stimulating hormone receptor Short name= FSH-R Alternative name(s): Follitropin receptor
<b>Expression Region</b>	18-687
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.