



Recombinant Human Adhesion G protein-coupled receptor L4 (ADGRL4)

Product Code	CSB-CF872541HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9HBW9
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>QNCTKTPCLPNAKCEIRNGIEACYCNMGFSGNGVTICEDDNECGNLTQSCGE NANCTNTE GSYYCMCVPGFRSSSNQDRFITNDGTVCIENVNANCHLDNVCIAANINKTLTKI RSIKEP VALLQEVYRNSVTDLSPTDIITYIEILAESSLLGYKNNTISAKDTLSNSTLTEFVK TVN NFVQRDTFVWVKLSVNHRRTHLTKLMHTVEQATLRISQSFQKTTEFDTNSTD IALKVFF FDSYNMKHIHPHMNDGDYINIFPKRKAAYDSNGNVAVAFVYYKSIGPLLSSS DNFLLKP QNYDNSEEEERVISSVISVSMSSNPPTLYELEKITFTLSHRKVTD RYRSLCAFW NYSPDT MNGSWSSEGCELTYSNETHTSCRCNHLTHFAILMSSGPSIGIKDYNILTRITQL GIIISL ICLAICIFTFWFFSEIQSTRTTIHKNLCCSLFLAELVFLVGINTNTNKLFC SIIAGLL HY FFLAFAWMCIEGIHLYLIVVGVYIYKNGFLHKNFYIFGYLSPAVVVGFSAALGYR YYGTT KVCWLSTENNF IWSFIGPACL IILVNLLAFGVIIYK VFRHTAGLKPEVSCFENIRS CARG ALALLFLLGTTWIFGVLHVHASVVTAYLFTVSNAFQGMFIFLFCVLSR KIQEE YYRLF KNPCCFGCLR</p>
Source	in vitro E.coli expression system
Target Names	ADGRL4
Protein Names	Recommended name: EGF, latrophilin and seven transmembrane domain-containing protein 1 Alternative name(s): EGF-TM7-latrophilin-related protein Short name= ETL protein
Expression Region	20-690
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.