



# Recombinant Mouse Lutropin-choriogonadotropic hormone receptor (Lhcgr)

<b>Product Code</b>	CSB-CF012911MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P30730
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Sequence</b>	<p>PELSGSRCPGPCDCAPDGALRCPGPRAGLARLSLTYPVKVIPSQAFRGLNEV VKIEISQ SDSLERIEANAFDNLNLSSEILIQNTKNLLYIEPGAFTNLPRLKYLSICNTGIRTL DVS KISSEFNFILEICDONLYITTIPGNAFQGMNNEITLKYGNGFEEVQSHAFNGTT LISL ELKENIYLEKMHSGETFQGATGPSILDVSSTKLQALPSHGLESIQTLIATSSYS LK TLPSR EKFTSLLVATLTYP SHCCAFRNLPKKEQNF SFSIFENFSKQCESTVREANN ETL YSAIFE ENELSGWDYDYDFCSPKTLQCTPEPDAFNPCEDIMGYAFLRVLIWLNILAIF G NLTVLF VLLTSRYKLTVPRFLMCNLSFADFCMGLYLLLIASVDSQTKGQYYNHAI DWQT GSGCSAA GFFT VFASELSVYTLTVITLERWHTITYAVQLDQKLRLRHAIPIMLGGWIF STLM ATLPL VGVSSYMKVSIKLPMDVESTLSQVYILSILLNNAFAFVVICACYVRIYFAV QNPEL TAPN KDTKIAKKMAILIFTDFTCMAPISFFAISA AFKVPLITVTNSKVLLVLFY PVNSCAN PFL YAVFTKAFQRDFLLLSRFGCKHRAELYRRKEFSACTFNSKNGFPRSSKPS Q AALKLSI VHCQQPTPPRVLIQ</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Lhcgr
<b>Protein Names</b>	Recommended name: Lutropin-choriogonadotropic hormone receptor Short name= LH/CG-R Alternative name(s): Luteinizing hormone receptor Short name= LSH-R
<b>Expression Region</b>	27-700
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.