



Recombinant Chicken Follicle-stimulating hormone receptor (FSHR)

Product Code	CSB-CF009021CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79763
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	QHHTCLCEGRIFICQEIKVVQLPRDIPTNATELRFVLTkmrvipkgAFTGLHDLE KIEIS QNDALeIIeGNVFSSLPKLHEIRIEKANKLMKIDQDAFQHLPsLRyLLISNTGLsFLP LPVV HKVHSFQKVLLDVQDNIHIRTIERNTFMGLSSESvILRLNkNGIQEIKDHAFNGT CLDEL NLSDNYNLEKLPEKVFQGAIGPVVLDISRTRISFLPSHGLeFIKkLRARSTYKLK KLpDV NKFRSLIEANFTYPSHCCAFTNRKTQNTeFYPICSMSPAKQDLGEQTGKRKHR RSAAEDY ISHYGTRFGPVENEFDYGLCNEVVDFVCSPKPDafNPCEdIMGYNVLRVLIWFI NILAIT GNTTVLIILISSQYKLTVPRLMCNLAFADLCIGIYLLFIASVDIQTKSRYYNYAID WQT GAGCNAAGFFTVFASeLSVYTLTVITLERWHTITYAMQLNrkVRLRHAVIIMVF GWMFAF TVALLPIFGISSYMKVSICLPMHIETPFsQAYVIFLLVLNVLAfVIICICYICIYFTVR N PNVISSNSDTKIakRMAILIFTDFLCMAPISFFAISASLRVPLITVSKSKILLVLFYPI N SCANPFLYAIFTKtFRRDFFILLsKFGCCeMQAQIYRTETSSSAHNfHTRNGHY PTASKN SDGTIYSLVPLNHLN
Source	in vitro E.coli expression system
Target Names	FSHR
Protein Names	Recommended name: Follicle-stimulating hormone receptor Short name= FSH-R Alternative name(s): Folitropin receptor
Expression Region	19-693
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.