



Recombinant Macropus eugenii Follicle-stimulating hormone receptor (FSHR)

Product Code	CSB-CF754488MPE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6YNB6
Product Type	Transmembrane Protein
Immunogen Species	Macropus eugenii (Tamar wallaby) (Notamacropus eugenii)
Sequence	CHHRICHCSDRVFCQESKVTEIPSDIPRNTVEMRFVLTKLQVIPRGAFSGFRDL EKIEI SQNDALEVIEADVFSNLPKLHEIRIEKANLVLIDPEAFWNLPNLRYSNTGIK HLP VYKIQSHQKVLLDIQDNINIRTIERNFSAGLSSESEILWLNKNGIQEIQNQAFNGT QLDE LNLSDNINLEDLPNGIFQGANGPVILDISRTIRHSLPSDGLKNLKKLKARSAYNF KTLPN LDKFAELVEANLTYPSHCCAFANWRRQAFELHPICNKSFVRQDSDDLAMDWG HRRRSVEE EYVSNLDKGFDTITDTELDYDLCNEVVDVACSPKPDFAFNPCEDIMGYNFLRVLIW FISLS ITGNIVVLVILITSQYKLTVPFLMCNLAFADLCIGIYLLIASVDIHTKSQYHNYAI DW QTGAGCDAAGFFTVFASLSVYTLTAITLERWHTITYAMQLDRKVRLRHAASIM LIGWIF AFSVALLPIFGVSSYMKVSICLPMDIDSPLSQFYVISLLVLNVLASVIICTCYTHIY FTV RNPNIISSTSDAKIAKRMAMLIFTDFLCMAPISFFAISASVKMPLITVSKSKILLVLF YP INSCANPFLYAVFTKTFRRDFFILLSKFGCCEMQAQIYRTETSSTVHSSHPRNG QCPLMA KNSGAVYKLVPLNHL
Source	in vitro E.coli expression system
Target Names	FSHR
Protein Names	Recommended name: Follicle-stimulating hormone receptor Short name= FSH-R Alternative name(s): Follicle-stimulating hormone receptor
Expression Region	18-694
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.