



Recombinant Pongo abelii Frizzled-6 (FZD6)

| | |
|--------------------------|--|
| Product Code | CSB-CF716182PYX |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | Q5RCN4 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii) |
| Sequence | <p>HSLFTCEPITVPRCMKMAYNMTFFPNLMGHYDQSIAAVEMEHFLPLANLECS NIETFLC KAFVPTCTEQIHVVPPCRKLCEKVYSDCKKLIDTFGIRWPEELEC DRLQYCD ETVPVTFD PHTEFLGPQKKTEQVQRDIGFWCPRHLKTS GGQGYKFLGIDQCAPP CPNMYF KSDELEFA KSFIGTVSIFCLCATLFTFLTLFLIDVRRFRYPERPIIYYSVCYSIVSLMYFIGFLLGD ST ACNKADEKLELGDTVVLGSQNKACTVLFMLLYFFTMAGTVWWVILTITWFLAA GRKWSCE AIEQKAVWFHAWAWGTPGFLTVMLLAMNKVEGDNISGVCFVGLYDL DASRYF VLLPLCLC VFGLSLLLAGIISLNHVRQVIQHDGRNQEK LKKFMIRIGVFSGLYLVPLVTLG CYVYE QVNRITWEITWVSDHCRQYHIPCPYQAKAKARPELALFMIKYLMTLIVGISAVF WVGSKK TCTEWAGFFKRNKRDPISERRVLQESCEFFLKHNSKVKHKKKHYKPS SHKL KVISKSM GTSTGATANHGTSVAITSHDYLGQETL TEIQTSPETSMREVKADGASTPRLR EQDCGEP ASPAASISRLSGEQVDRKGQAGSVSE SARSEGRISPKSDITDTGLAQSNNLQV PSSSEPS SLKGSTSLVHPVSGVRKEQGGGCHSDT</p> |
| Source | in vitro E.coli expression system |
| Target Names | FZD6 |
| Protein Names | Recommended name: Frizzled-6 Short name= Fz-6 |
| Expression Region | 20-707 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |