



# Recombinant Human Frizzled-6 (FZD6)

<b>Product Code</b>	CSB-CF009109HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O60353
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sequence</b>	<p>HSLFTCEPITVPRCMKMAYNMTFFPNLMGHYDQSIAAVEMEHFLPLANLECS  NIETFLC  KAFVPTCIEQIHVVPPCRKLCEKVYSDCKKLIDTFGIRWPEELEC  DRLQYCD  ETVPVTFD  PHTEFLGPQKKTEQVQRDIGFWCPRHLKTSGGQGYKFLGIDQC  APPCPNMYF  KSDELEFA  KSFIGTVSIFCLCATLFTFLTLFLIDVRRFRYP  ERPIIYYSVCYSIVSLMYFIGFLLGD  ST  ACNKADEKLELGDTVVLGSQNKACTVLFMLLYFFTMAGT  VWWWVILTITWFLAA  GRKWSCE  AIEQKAVWFHAWAWGTPGFLTVMLLAMNKVEGD  NISGVCFVGLYDL  DASRYFVLLPLCLC  VFGLSLLLAGIISLNHVRQVIQHDGRNQEK  LKKFMIRIGVFSGLYLVPLVTLG  CYVYE  QVNRITWEITWVSDHCRQYHIPCPYQAKAKARPE  LALFMIKYLMTLIVGISAVF  WVGSKK  TCTEWAGFFKRNKRDP  ISESRRVLQESCEFFLKHNSKV  KHKKKHYKPSSHKL  KVISKSM  GTSTGATANHGTS  AVAITSHDYLGQETL  TEIQTSPETSMRE  VKADGASTPRLR  EQDCGEP  ASPAASISRLS  GEQVDGKGQAGSV  SEARSSEGRISPK  SDITDTGLAQS  NNLQV  PSSSEPS  SLKGSTSLVHPV  SGVRKEQGGG  CHSDT</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	FZD6
<b>Protein Names</b>	Recommended name: Frizzled-6 Short name= Fz-6 Short name= hFz6
<b>Expression Region</b>	19-706
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.