



Recombinant Mouse Frizzled-6 (Fzd6)

Product Code	CSB-CF717196MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q61089
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p>HSLFTCEPITVPRCMKMTYNMTFFPNLMGHYDQGIAAVEMGHFLHLANLECS NIEMFLC QAFIPTCTEQIHVVLPCKRLCEKIVSDCKKLMDFGIRWPEELECNRPLPHCDDT VPVTSH PHTELSGPQKKSDQVPRDIGFWCPKHLRTSGDQGYRFLGIEQCAPPCPNMYF KSDEL DFA KSFIVSIFCLCATLFTFLTLIDVRRFRYPERPPIIYYSVCYSIVSLMYFVGFLLGN ST ACNKADEKLELGDTVVLGSKNKACSVVMFLYFFTMAGTVWWVILTITWFLAA GRKWSCE AIEQKAVWFHAWAWGAPGFLTVMLLAMNKVEGDNISGVCFVGLYDL DASRYF VLLPLCLC VFGLSLLLAGIISLNHVRQVIQHDGRNQEKLLKFMIRIGVFSGLYLVPLVTLG CYVYE LVNRTWEMTWFS DHCHQYRIPCPYQANPKARPELALFMIKYLMTLIVGISAVF WVGSKK TCTEWAGFFKRNRKRDPISERRVLQESCEFFLKHNSKVKKHKKHKGAPGPHR LKVISKSM GTSTGATTNHGTSAMAIADHDYLGQETSTEVHTSPEASVKEGRADRANTPSA KDRDCGES AGPSSKLSGNRNGRESRAGGLKERSNGSEGAPSEGRVSPKSSVPETGLIDCS TSQAASSP EPTSLKGSTSLPVHSASRARKEQGAGSHSDA</p>
Source	in vitro E.coli expression system
Target Names	Fzd6
Protein Names	Recommended name: Frizzled-6 Short name= Fz-6 Short name= mFz6
Expression Region	19-709
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.