



# Recombinant Bovine CD97 antigen (CD97)

<b>Product Code</b>	CSB-CF823366BO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q8SQA4
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Sequence</b>	<p>QNSKACALPCPPNSSCVNGTACRCAPGFISFSGEIFTDPLESCDDINECGPPS  PVDCGSS  ADCQNTEGGYYCTCSPGYEPVSGAMIFRNESENTCRDVDECSSGQHQCCHNS  TVCFNTVGS  YTCHCREGWEPKHGLKNKQKDTICKEISFPAWTAPPGIKSRSLSAFFERVQKM  SRDFKPA  MAKXSMQDLVGSVDDLLKNSGDLES LDQSSKHVTVTHLLSGLEQILRTLAKAM  PKGSFTY  RSLDNTELSLVVQEQGKGNVTVGQSHARMLLDWAVAAAAEESGPTVVGILSS  QNMKKLLA  NASLKL DSEKLYKSPVRGAKVTLLSAVSSVFLSNTNTEKLD SNV SFAFALH  EQPELK  PRQELICAFWKKDSNGNGSWATTGCWKMGNGSITCQC SHLSSFAILMAH  YDVEDPKLA  LITKVG LALSLACLLLCILTFLLVRPIQGSRTTVHLHLCICL FVGS AIFLAGIENEG  GEV  GTRCRLVAVLLHYCF LA AFCWMSLEGVELYFLVVRV FQGQGM RKLWLCLIGY  GVPLIIVG  ISAGAYSKGYGREKFCWLNFE GGF LWSFVGPVTFIVLGN AIIFVITVWKL TQKF  SEINPD  IKKLLKARVLTITAI AQLFVLGCTWV FGLLLFN PESWVLSYIF SILNCLQGFFLV L  YCL LNKKVREEYR KWACMVAGNKYSEFATTTSGSGSSH NQTQALRPSESGM</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	CD97
<b>Protein Names</b>	Recommended name: CD97 antigen Alternative name(s): CD_antigen= CD97
<b>Expression Region</b>	27-734
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.