



Recombinant Mouse EGF, latrophilin seven transmembrane domain-containing protein 1 (EltD1)

Product Code	CSB-CF835658MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q923X1
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p>QNCSTTCLPNAKCEVHNGVEACFCSQGYSGNGVTICEDIDECSESSVCGDH AVCENVNG GFSCFCREGYQTATGKSQFTPNDGSYCQDIDECSESSVCGDHAVCENVNGG FSCFCREGY QTATGKSQFTPNDGSYQCESMNSNCHLEHACIAANINKTLKRIGPITEQTLLQ EIYRNS EAELSLMDIVTYIEILTESSLLGHPNSTTSYKDAHFNSTLTEFGETINNFVERST HKMW DQLPTNHRRLHLTKLMHTAELVTLQIAQNTQKNSQFDMNSTDLALKVFAFDST HMKHAHP HMNVDDGGYVKISPRRKAAGTTGNVVVAFLCYKSIGPLLSSSDNFLLDTQNDN SEGKEKV ISSVISASISSNPPTLYELEKITFTLSHVKLSDKHRTQCAFWNYSVDAMNNGSW STEGCE LTHSNDTHTSCRCSHLTHFAILMSSTSSIGIKDYNILTRITQLGIIISLICLAICITF W FFSEIQSTRTTIHKNLCCSLFLAELVFLIGININTNKLVCIIAGLLHYFFLAFAW MCI EGIHLYLIVVGVYINKGFLHKNFYIFGYLSPA VVVGFSASLGYRYYGTTKVCWLS TENNF IWSFIGPACLIILVNLLAFGVIIYKVRHTAGLKPEVSCYENIRSCARGALALLFLL GTT WIFGVLHVHVASVVTAYLFTVSNAFQGMFIFLFLCVLSRKIQEEYYRLFKNVPC CFGCLR</p>
Source	in vitro E.coli expression system
Target Names	Adgrl4
Protein Names	Recommended name: EGF, latrophilin seven transmembrane domain-containing protein 1
Expression Region	20-739
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.



Protein Length

Full Length of Mature Protein

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.