



# Recombinant Mouse Relaxin receptor 1 (Rxfp1)

<b>Product Code</b>	CSB-CF764847MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q6R6I7
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Sequence</b>	<p>MTSGPFFFCIFIIGKYFTLGSAQDVSCPLGSFPCGNMSRCLPQLLHCNGVDDC  GNRADED  HCGDNNGWSLQLDKYFANYKLASTNSFEAETSECLVGSVPMHCLCRDLELD  CDEANLRA  VPSVSSNVTVMQLQRNFIRTLPPNGFRKYHELQKLCLQNNRIHSVSVSAFRGL  RSLTKLY  LSHNRITFLKPGVFEDLHRLEWLIIEDNHLRSRISPLTFYGLNSLILLVLMNNALTR  LPDK  PLCQHMPRLHWLDFEGNRIHNLRLNLTIFSCNNLTVLVMRKNKINYLNEHAFTHL  QKLDEL  DLGSNKIENLPPNIFKDLKELSQLNISYNPIQKIEVNQFDCLAKLKSLSLEGIEISN  IQQ  RMFRPLINLSHIYFKKFQYCGYAPHVRSCKPNTDGISSLENLLASIIQRVFWVVV  SAITC  FGNIFVICMRPYIRSENKLHAMSIIISLCCADCLMGVYLFVIGAFDLKFRGEYNKH  AQPWM  ESVHCQFMGSLAILSTEVSLLLLTFLTLEKYICIVYFPRCLRPRKCRITVLIWII  GF  IVAFAPLGNKEFFKNYYGTNGVCFPLHSEDTGSTGAQIYSVIFLGINLVAFIIIVF  SYG  SMFYSVHQSSVTVTEIQKQVKKEVVLAKRFFIVFTDALCWIPFILKFLSLLQVEI  PDS  ITSWVIFILPINSALNPIIYTLTRPFKEMIHQLWHNYRQRRSVDRKETQKAYAP  SFIW VEMWPLQEMSSGFMKPGAFTDPCDLSLVSQSSRLNSYS</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Rxfp1
<b>Protein Names</b>	Recommended name: Relaxin receptor 1 Alternative name(s): Leucine-rich repeat-containing G-protein coupled receptor 7 Relaxin family peptide receptor 1
<b>Expression Region</b>	1-758
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.