



Recombinant Human Probable G-protein coupled receptor 128 (GPR128)

Product Code	CSB-CF853445HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96K78
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>LGIWRIVIRIQRGKSTSSSSTPTEFCRNGGTWENGRICICTEEWKGLRCTIANFC ENSTYM GFTFARIPVGRYGPSLQTCGKDTPNAGNPMVRLCSLSLYGEIELQKVTIGNC NENLETL EKQVKDVTAPLNNISSEVQILTSDANKLTAENITSATR VVGQIFNTRNASPEAK KVAIV TVSQLLDASEDAFQRVAATANDDALTTLIEQMETYSLSLGNQSVVEPNIAIQSA NFSSSEN AVGPSNVRFSVQKGASSSLVSSSTFIHTNVDGLNPDAQTELQVLLNMTKNYTK TCGFVVY QNDKLFQSKTFTAKSDFSQKIISSKTDENEQDQSASVDMVFSPKYNQKEFQLY SYACVYW NLSAKDWDTYGCQKDKGTDGFLRCRCNHTTNFAVLMTFKKDYQYPKSLDILS NVGCALSV TGLALTVIFQIVTRKVRKTSVTWVLVNLCSMLIFNLLFVFGIENSKNLQTSDDG INNI DFDNDIPRTDTINIPNPMCTAIAALLHYFLLVTFTWNALSAAQLYLLIRTMKPL PRHF ILFISLIGWGVPAIVVAITVGVVIYSQNGNPNQWELDYRQEKICWLAIPEPNGVIKS PLLW SFIVPVTIILISNVVMFITISIKVLWKNNQNLSTKKVSSMKKIVSTLSVAVVFGITW IL AYLMLVNDD SIRIVFSYIFCLFNTTQGLQIFILYTVR TKVFQSEASKVLM LLSIG RRKS LPSVTRPRLRVKMYNFLRSLPTLHERFR LLETSPSTEEITLSESDNAKESI</p>
Source	in vitro E.coli expression system
Target Names	ADGRG7
Protein Names	Recommended name: Probable G-protein coupled receptor 128
Expression Region	27-797
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.