



Recombinant Rat Probable G-protein coupled receptor 156 (Gpr156)

Product Code	CSB-CF817004RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8K451
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>MEPEINCSEFCDSFPGQELDRRPLHDLCKTTITDSQHGSADISPLSPALLGVIW TFLSCG LLLVLFFLAFTIRCRKNRIVKMSSPNLNIVTLLGSCCLTYSSAYLFGIQDALVGSSV EALI QTRLSELLCIGTTLVFGPILGKSWRLYKVFTQRVDPKRVIIKDLQLLGLVAALVVA DVILL VTWVLTDPICLQILGVSMKVTGRDVSCSLTNTTHFCASRYSDVWIALVLGCKG LLLLYGA YLAGLTHVSSPPVNQSLTIMVGVNLLLLTAGLLFVVTRYLHWPNLVFGLTSG GIFVCT TTVNCCVFLPQLRQRKAFEGENQTIHMAKYFSTPSKTFRSKFDEDQSCHLR DEXSCMER LLTEKNAVIESLQEQVSNAKEKLVKLMSAECALDSPEWAVPAAASAGGPAECX ATSEKES GAAAEDSLPASAASQHMQGPGASRRDXSPSPDQKYDMPLKQFCDHLDMGC SQKPKAEQSE GPERGNQEPMAPGQSLMTDGVACEPHRPRQNSEVLPERLPRVSSVVREKLQ EVLQELDLG TEAPLSPLPCPQQPWKSNTSGSPQKLSKLGFSYVRRRRRAAQRARSRIP GSVGLKMG HQANNTVSGSQNGLIVQNRDSPRLDHHNARSKEPRSSSVKPSISPAPHQRRG SLEGSKQC ETEPQEARGYSVAFPRQPSASAPAQSSTAPCLSSXPALPRQRQPLPLLSPGC PSLSSGCY NLDSESSSSDEFFCRCHRPYCEICFQSSLDSDNDSDTSDSDLEQTSGLASWEK LLARSKPV VNFKDDLKPTLV</p>
Source	in vitro E.coli expression system
Target Names	Gpr156
Protein Names	Recommended name: Probable G-protein coupled receptor 156 Alternative name(s): GABAB-related G-protein coupled receptor
Expression Region	1-792
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.



Protein Length

full length protein

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.