



Recombinant Mouse Probable G-protein coupled receptor 156 (Gpr156)

Product Code	CSB-CF761328MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6PCP7
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p>MEPEINCSEFCDSFPGQELDRRPLHDLCKTTITESQHSSTAASPLSPALLGIMW TFLSCG LLLVLFFLAFTIRCRKNRIVKMSSPNLNVVTLGSCCLTYISAYLFGIQDALEGSSV EALI QTRLSELLCIGTSLVFGPILGKSWRLYKVFTQRVPDKRVIKDLQLLGLVAALVVA DVILL VTWVLTDPICLQMLGVSMKVTGRDVSCSLTNTHFCASRYSDVWIALVLGCK GLLLLYGA YLAGLTHVSSPPVNQSLTIMVGVNLLLLTAGLLFVVTRYLHWPNLVFGLTSG GIFVCT TTVNCCVFIPQLKQWKAFEGENQTMRHMAKYFSTPSKSFHSQFDEDPSCHLR DEKSCMER LLTEKNAVIESLQEQVSNAKEKLVKLMSAECTYDSPEWAVPDAASARGLALPG PSECPAV SENEGAAARDSLHVPAACQHVQGPASRRDTSPSPAQQDNMPLKQYCDHL DTGCNQKPK AEQSEGPBERGDQEPMAPSQRLMADGVACEPHKPRQSPEGLPKKLPGVSSV REKLQEVLQ ELDLGSEAPLSPLPCPQQLWKSTTSRSPQKLSKLGFSFYVRRRRRAAQR RSHIPGSV GLNVGHQANSTVSSSQSGLIVQNRDSPRLDHHNARSKVPRSSSVKPSPLSEP RRKQGTLE GSKQCETEPQEAGGACNVAFFPCQSSASVQAQSPAAPCLPSSPALPRQRQPR PRLSPGCPS LSSGCYNLDSESSSSDEFFCRCHRPYCEICFQSSLDSNDSDTSDSDLEQASGL ASWEKLW ARSKPVVNFKDDLKPTLV</p>
Source	in vitro E.coli expression system
Target Names	Gpr156
Protein Names	Recommended name: Probable G-protein coupled receptor 156 Alternative name(s): GABAB-related G-protein coupled receptor
Expression Region	1-798
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.



Protein Length

full length protein

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.