



Recombinant Human EGF-like module-containing mucin-like hormone receptor-like 2 (EMR2)

Product Code	CSB-CF887050HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9UHX3
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>QDSRGCARWCPQDSSCVNATACRCNPGFSSFSEIITTPMETCDDINECATLSK VSCGKFS DCWNTGEGSYDCVCSPGYEPVSGAKTFKNESENTCQDVDECQQNPRLCKSYG TCVNTLGSY TCQCLPGFKLPEDPKLCTDVNECTSGQNPCHSSTHCLNNVGSYQCRCRPG WQPIPGSPN GPNNTVCEDVDECSSGQHQCDSSTVCFNTVGSYSCRCRPGWKPRHGIPNN QKDTVCEDMT FSTWTPPPGVHSQTLRFFDKVQDLGRDYKPGLANNTIQSILQALDELLEAPG DLETLPR LQQHCVASHLLDGLLEDVLRGLSKNLSNGLLNFSYPAGTELSLEVQKQVDRSVT LRQNQAV MQLDWNQAQKSGDPGPSVVGLVSIPGMGKLLAEAPLVLEPEKQMLLHETHQ GLLQDGSP LLSDVISAFLSNNDTQNLSSPVTFTFSHRVIPRQKVLVCFWEHGQNGCGHWA TTGCSTI GTRDTSTICRCTHLSSFAVLMAYHDVQEEDPVLTVITYMGLSVSLLCLLLAALTF LLCKA IQNTSTSLHLQLSLCLFLAHLFLVAIDQTGHKVLCSIIAGTLHYLYLATLTWMLL EALY LFLTARNLTVVNYSSINRFMKKLMFPVGYGPAVTVVAISAASRPHLYGTPSRC WLQPEKG FIWGFLGPVCAIFSVNLVFLVTLWILKNRSSLNSEVSTLRNTRMLAFKATAQL FILGC TWCLGILQVGPAARVMAYLFTIINSLQGVFIFLVYCLLSQQVREQYGKWSKIR KLKTES EMHTLSSSAKADTSKPSTVN</p>
Source	in vitro E.coli expression system
Target Names	ADGRE2
Protein Names	Recommended name: EGF-like module-containing mucin-like hormone receptor-like 2 Alternative name(s): EGF-like module receptor 2 CD_antigen=CD312
Expression Region	24-823
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.



Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>