



Recombinant Danio rerio Leukocyte cell-derived chemotaxin 1 (lect1)

Product Code	CSB-CF350134DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P58239
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	SATRMRRQTSAGVNRQPARRRNSTASARDERPTGPEYNPENPYHQNQGSE GTMVFDPM LD HRGICCTECHRSYTHCERVCEPLGGYWPWPYNYHGCRPVCR LIMP CRWWA ARVLGLV
Source	in vitro E.coli expression system
Target Names	cnmd
Protein Names	Recommended name: Leukocyte cell-derived chemotaxin 1 Cleaved into the following 2 chains: 1. Chondrosurfactant protein Short name= 2. CH-SP 3. Chondromodulin-1 Alternative name(s): Chondromodulin-I Short name= ChM-I
Expression Region	170-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.