



Recombinant Rabbit ATP synthase protein 8 (MT-ATP8)

Product Code	CSB-CF015071RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O79431
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	MPQLDTSTWFTTIVAMILSLFILMQLKFHKYTYPMNPVLKALESTSFPCPWETK WTKIYSPLSLPQH
Source	in vitro E.coli expression system
Target Names	MT-ATP8
Protein Names	Recommended name: ATP synthase protein 8 Alternative name(s): A6L Chargerin II F-ATPase subunit 8
Expression Region	1-67
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.