



Recombinant *Cricetulus griseus* ATP synthase protein 8 (MT-ATP8)

| | |
|--------------------------|--|
| Product Code | CSB-CF015071DXU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P14414 |
| Product Type | Transmembrane Protein |
| Immunogen Species | <i>Cricetulus griseus</i> (Chinese hamster) (<i>Cricetulus barabensis griseus</i>) |
| Sequence | MPQLDTSTWFTTVLASSITLFILMQLKISFHDLHKKPSNKYLKLFKPTNPWEQK WTKIYSPLSLPQP |
| Source | in vitro E.coli expression system |
| Target Names | MT-ATP8 |
| Protein Names | Recommended name: ATP synthase protein 8 Alternative name(s): A6L F-ATPase subunit 8 |
| Expression Region | 1-67 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |