



Recombinant *Artibeus jamaicensis* ATP synthase protein 8 (MT-ATP8)

Product Code	CSB-CF015071DOQ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O99598
Product Type	Transmembrane Protein
Immunogen Species	<i>Artibeus jamaicensis</i> (Jamaican fruit-eating bat)
Sequence	MPQLDTSTWFITILATILTLFIIMQLKISTYYYHSNPEPKTTKMTKSLIPWEIKWTKI YSPLSLPLR
Source	in vitro E.coli expression system
Target Names	MT-ATP8
Protein Names	Recommended name: ATP synthase protein 8 Alternative name(s): A6L F-ATPase subunit 8
Expression Region	1-67
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.