



Recombinant Dog Phospholemman (FXVD1)

Product Code	CSB-CF009089DO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P56513
Product Type	Transmembrane Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Sequence	EAPQEHDPTFDYQSLRIGGLIIAGILFILGILIVLSRRRCRCKFNQQQRTGEPDEE EGTFRSSIRRLSTRRR
Source	in vitro E.coli expression system
Target Names	FXVD1
Protein Names	Recommended name: Phospholemman Alternative name(s): FXVD domain-containing ion transport regulator 1
Expression Region	21-92
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.