



Recombinant Human Putative FXYP domain-containing ion transport regulator 8 (FXYP6P3)

Product Code	CSB-CF009098HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P58550
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	SAAKEKEIDPFHYNQTLRIGGLVFDVVLFLVPSCHLLSHRCKCSFNQKPQDP GDKEAQVENFITANAKEPQKAKN
Source	in vitro E.coli expression system
Target Names	FXYP6P3
Protein Names	Recommended name: Putative FXYP domain-containing ion transport regulator 8
Expression Region	19-94
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.