



Recombinant Human TYRO protein tyrosine kinase-binding protein (TYROBP)

Product Code	CSB-CF025398HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43914
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	LRPVQAQAQSDCSCSTVSPGVLAGIVMGDLVLTVLIALAVYFLGRLVPRGRGA AEAATRKRITETESPYQELQGQRSDVYSDLNTQRPYYK
Source	in vitro E.coli expression system
Target Names	TYROBP
Protein Names	Recommended name: TYRO protein tyrosine kinase-binding protein Alternative name(s): DNAX-activation protein 12 Killer-activating receptor-associated protein Short name= KAR-associated protein
Expression Region	22-113
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.