



Recombinant Dog Protein transport protein Sec61 subunit beta (SEC61B)

Product Code	CSB-CF020958DO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P60467
Product Type	Transmembrane Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Sequence	PGPTPSGTVNGSSGRSPSKAVAARAAGSTVRQRKNASCGRTRSAGRTTSAGT GGMWRFYTEDSPGLKVGPPVLMVMSLLFIASVFMLHIWKGKYTRS
Source	in vitro E.coli expression system
Target Names	SEC61B
Protein Names	Recommended name: Protein transport protein Sec61 subunit beta
Expression Region	2-96
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.