



Recombinant Mouse BET1-like protein (Bet1l)

Product Code	CSB-CF002671MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O35153
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MADWTRAQSSGAVEDILDRENKRMADSLASKVTRLKSLALDIDRD TEDQNRYL DGMDSDFTSVTGLLTG SVKRFSTMARSGRDNRKLLCGMAVVLIVAFFILSYLLS RTRT
Source	in vitro E.coli expression system
Target Names	Bet1l
Protein Names	Recommended name: BET1-like protein Alternative name(s): Golgi SNARE with a size of 15 kDa Short name= GOS-15 Short name= GS15 Vesicle transport protein GOS15
Expression Region	1-111
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.