



Recombinant Human Cell cycle control protein 50C (TMEM30C)

Product Code	CSB-CF023827HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A0ZSE6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MEERAQHCLSRLLDNSALKQQELPIHRLYFTARRVLFVFFATGIFCLCMGIILILS ARSTQEIEINYTRICANCAKLRENASNFDKECTCSIPFYLSGKMMVGEIQETRLT LH
Source	in vitro E.coli expression system
Target Names	TMEM30CP
Protein Names	Recommended name: Cell cycle control protein 50C Alternative name(s): Transmembrane protein 30C
Expression Region	1-113
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.