



Recombinant Bovine Vesicle-associated membrane protein 2 (VAMP2)

Product Code	CSB-CF025781BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P63026
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	SATAATAPPAAPAGEGGPPAPPNLTSNRRLQQTQAQVDEVVDIMRVNVDKV LERDQKLSELDDRADALQAGASQFETSAAKLKRKYWWKNLKMMLGVICAILII IIVYFSS
Source	in vitro E.coli expression system
Target Names	VAMP2
Protein Names	Recommended name: Vesicle-associated membrane protein 2 Short name= VAMP-2 Alternative name(s): Synaptobrevin-2
Expression Region	2-116
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.